



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,
Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Ice Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC040101
Material	Bandage Composition: a. Blue Elastic Bandage: 65% polyester, 15% spandex and 20% cotton b. Non Woven Cohesive Bandage: Non Woven Fabric, Latex Glue or Latex Free Synthetic Glue c. Cooling Ingredients: distilled water, medical alcohol, mint and mixed borneol
Color	Blue/Colored Cohesive Bandage Colors
Size	Width: 7.5cm, 8cm, 10cm Length (Stretched): 3.0m, 3.6m
Package	1 Roll/Pack
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Ice Bandage Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package	Carton Size
JC0401014	Ice Bandage, instant cooling relief, 7.5/8.0/10cm x 3.6m, 1 roll/pack.	7.5cm x 3.6m	1 Roll/Pack	60 Packs	43*20*30CM
JC0401015		8cm x 3.0m	1 Roll/Pack	60 Packs	43*20*30CM
JC0401016		10cm x 3.6m	1 Roll/Pack	48 Packs	43*20*30CM